Biochemistry

Acid Phosphatases (Total & Prostatic)(Manual)	01
Alanine Aminotransferase (ALT, GPT)(Automated/Manual)	01
Aspartate Aminotransferase (AST, GOT)(Automated/Manual)	01
Aminotransferase (AST & ALT)(Manual)	01
Albumin (Manual)	01
Alkaline Phosphatases (Automated/Manual)	01
Alkaline Phosphatases (Manual)	01
α-Amylase (Manual)	01
Bilirubin (Total; Automated/Manual)	01
Bilirubin (Direct; Automated/Manual)	01
Bilirubin (Total & direct; Manual)	01
Calcium (Manual)	01
Calcium (Automated/Manual)	01
Cholesterol (Automated/Manual)	01
Copper (Automated/Manual)	01
Creatinine (Automated/Manual)	01
Creatinine (Manual)	01
Creatine Phosphokinase (CK)(Automated/Manual)	01
Glucose (Automated/Manual)	01
Glucose 6-Phosphate Dehydrogenase, G6PD (Colorimetric, Manual)	01
Hemoglobin A1c (Turbidimetric, Automated)	01
High Density Lipoprotein (HDL)- Cholesterol (Direct, Automated & Manual)	01
High Density Lipoprotein (HDL)- Cholesterol (Precipitation, Automated/Manual)	01
Low Density Lipoprotein (LDL)- Cholesterol (Direct)	01
Iron (Automated/Manual)	01
Iron Binding Capacity,TIBC (Direct, Automated/Manual)	01
Iron Binding Capacity (Pricipitation, Automated/Manual)	01
Lactate Dehydrogenase (LDH)(Automated/Manual)	01
Magnesium (Automated/Manual)	01
Phosphorus (Automated/Manual)	01
Phosphorus (Manual)	01
Total Protein (Automated/Manual)	01
Triglycerides (Automated/Manual)	01
Urea (UV, Automated/Manual)	01
Urea (Berthelot, Automated/Manual)	01
Urea (BUN; Manual)	01
Uric Acid (Automated/Manual)	01
Uric Acid (Manual)	01
Urine & CSF Protein(Micro-Turbidimetric; Automated)	01
Urine and CSF proteins (Colorimetric; Automated/Manual)	01
Urinary Citrate (Automated/Manual)	01
Urinary Citrate (Manual)	01
Urinary Oxalate (Automated/Manual)	01
Urinary Calculi (Manual)	01
Urinary Cystine (Manual)	01
Zinc (Automated/Manual)	01

Hematology

1. Diluents	01
2. Stromatolyser	01
3. Enzymatic Cleanser	01

IN THE NAME OF GOD

Manufacturing accurate and reliable biochemical diagnostic kits in Iran, has always been one of our big desire in the past. It became more and more necessary during a period in which access to foreign made kits proved to be most difficult and in some cases rather impossible.

Darman Faraz kave diagnostic establishment was the result of the above idea and we do hope that we will soon be in such position not only to fulfill the national needs for biochemical diagnostics but to export to the neighbor countries providing government permission is obtained. We would like to express our sincere gratitude to our most respectful and knowledgable colleagues all over the country for their valuable advise, encouragements and criticisms, each of which has act as stimulation for us to pay more and more attention to our quality control department.

It is a greate honor for Darman Faraz kave to be able to announce that its quality control section is one of the best of its kind, comparable to any quality control department regardless its geographical position. Once again we would like to express our thankfulness to all governmental sections involved in diagnostics, to our friends and colleagues all over the country for their encouragements, advice and positive criticisms and we do look forwards to see Darman Faraz kave is presentable internationally.

Darman Faraz kave

### Acid Phosphatases (Total & Prostatic)

Method: Colorimetric, Manual

Principle:

The enzyme is permitted to act on p-nitrophenyl phosphate at pH 4.9, in both the presence and absence of 0.04M tartrate ion, using a 30-minute reaction time. The liberated p-nitrophenol is measured colorimetrically at 410nm after conversion to the quinoid form by sodium hydroxide. The intensity of the color is directly proportional to the concentration of the enzyme in the sample. The assay tube without tartrate gives a measure of all forms of the enzyme (platelet, red blood cell, prostatic etc). The tartrate inhibits the activity of the prostatic form and therefore the tube containing the tartrate measures the activity of the non-prostatic forms. The level of the prostatic enzyme is then obtained by difference.

Specimens:

Serum is the sample of choice. Plasma should not be used because many anticoagulants react with acid phosphatases.

Hemolyzed samples should not be used.

Contents:

Reagent 1: Buffer Solution (pH 4.9)	1 x 50ml
Reagent 2: PNPP substrate vials	10vials
Reagent 3: Tartrate reagent	1 x 6ml
Reagent 4: NaOH (5M)	1 x 5ml
Reagent 5: ACP standard	1 x 5ml
Number of Manual Assays:	50
Catalogue Number:	DFK-115

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### Alanine Aminotransferase (ALT, GPT)

Method: Enzymatic, UV Kinetic, Manual/Automated

Principle:

ALT in the sample catalyzes the transfer of an amino group from L-alanine to oxoglutarate. The pyruvate thus formed then reacts with nicotinamide-adenine dinucleotide (NADH) oxidizing it to NAD+. This reaction is catalyzed by the enzyme lactate dehydrogenase (LDH). The rate of decrease in the concentration of NADH is directly proportional to the activity of ALT and is monitored colorimetrically at 340nm

Specimens:

Serum and plasma give similar results although it is preferable to use serum. Slight hemolysis does not affect the result.

Contents:

Reagent 1: Enzyme Reagent	2 x 40ml
Reagent 2: Substrate Reagent	1 x 20ml
Number of Assays: Manual	100
Automated Dependent on Instrument type	
Catalogue Number:	DFK-128

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### Aspartate Aminotransferase (AST, GOT)

Method: Enzymatic, UV Kinetic, Manual/Automated

Principle:

AST in the sample catalyzes the transfer of an amino group from L-aspartate to  $\alpha$ -oxoglutarate. The oxaloacetate thus formed then reacts with nicotinamide-adenine dinucleotide (NADH) oxidizing it to NAD<sup>+</sup>. This reaction is catalyzed by the enzyme malate dehydrogenase (MDH). The rate of decrease in the concentration of NADH is directly proportional to the activity of AST and is monitored colorimetrically at 340nm. The enzyme lactate dehydrogenase (LDH) is incorporated in to the reagents in order to convert endogenous pyruvate in the sample to lactate during the lag phase of the analysis.

Specimens:

Serum and plasma give similar results although it is preferable to use serum. Hemolyzed samples should be avoided.

Contents:

Reagent 1: Enzyme Reagent	2 x 40ml
Reagent 2: Substrate Reagent	1 x 20ml
Number of Assays: Manual	100
Automated Dependent on Instrument type	
Catalogue Number:	DFK-126

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۲ × ۴۰ ml  
۲۰ × ۱ ml

### Aminotransferases

Aspartate Aminotransferase (AST, GOT)

Alanine Aminotransferase (ALT, GPT)

Method: Colorimetric, Manual

Principle:

AST in the sample catalyzes the transfer of an amino group from L-aspartic acid to  $\alpha$ -oxoglutaric acid. The oxalacetic acid thus formed can be assayed colorimetrically (505nm) by coupling with the chromogen 2,4-dinitrophenylhydrazine in the presence of sodium hydroxide. The intensity of the resultant colored complex is directly proportional to the activity of AST. Similarly, ALT in the sample catalyzes the transfer of an amino group from L-alanine to  $\alpha$ -oxoglutaric acid. The pyruvic acid thus formed is the assayed colorimetrically (505nm) by coupling with the chromogen 2,4-dinitrophenylhydrazine in the presence of sodium hydroxide. The intensity of the resultant colored complex is directly proportional to the activity of ALT.

Specimens:

Serum and plasma give similar results although it is preferable to use serum. Slight hemolysis does not affect the determination of ALT. However, hemolyzed sample should not be used for the determination of AST.

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Method: Colorimetric, Kinetic, Manual/Automated

Principle:

Alkaline phosphatase in the sample is permitted to act on p-nitrophenyl phosphate (PNPP) at pH 10.4 to produce p-nitrophenol (PNP) and inorganic phosphate. The yellow color of PNP is measured calorimetrically at 405nm. The rate at which PNP is produced is directly proportional to the activity of alkaline phosphatase.

Specimens:

Serum is the sample of choice. However, heparinized plasma can also be used. Many other anticoagulants inhibit alkaline phosphatase.

Contents:

Reagent 1: Buffer Solution (pH 10.4)	1 x 110ml
Reagent 2: PNPP Substrate	10vials
Number of Assays: Manual	100
Automated Dependent on the instrument type	
Catalogue Number:	DFK-129

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Method: Colorimetric, Manual

Principle:

Alkaline phosphatase in the sample is permitted to act on p-nitrophenyl phosphate (PNPP) at pH 10.4 to produce p-nitrophenol (PNP) and inorganic phosphate. The yellow color of PNP is measured calorimetrically at 410nm and is directly proportional to the activity of alkaline phosphatase.

Specimens:

Serum is the sample of choice. However, heparinized plasma can also be used. Many other anticoagulants inhibit alkaline phosphatase.

Contents:

Reagent 1: Buffer Solution (pH 10.4)	1 x 100ml
Reagent 3: NaOH (3M)	1 x 10ml
Reagent 2: PNPP Substrate	10 vials
Reagent 4: Stock Standard	1 x 5ml
Number of Manual Assays:	100
Catalogue Number:	DFK-112

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## Albumin

Method: Colorimetric, End-point, Manual

Principle:

Albumin is analyzed in the sample directly by the method of dye-binding. Bromocresol green binds readily to albumin and forms a green complex. The intensity of the color is directly proportional to the concentration of albumin and is measured colorimetrically at 630nm.

Specimens:

Serum is the sample of choice as plasma contains fibrinogen which, being a protein can interfere with the results.

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Method: Colorimetric, Manual

Principle:

The method is based on the ability of amylase to catalyze the hydrolysis of starch to simple sugars. The starch is incubated with serum and the rate of its disappearance (hydrolysis) is detected using iodine as an indicator.

Specimens:

Serum is the sample of choice. Hemolyzed sample should be avoided.

Contents:

Reagent 1: Substrate solution	2 x 50ml
Reagent 2: Stock iodine solution	1 x 10ml
Number of Manual Assays:	100
Catalogue Number:	DFK-113

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### Bilirubin (Total Automated)

Method: colorimetric , End Point , Automated and manual.

Principle: Congugated and unconjugated bilirubin (Total bilirubin) react with diazotized sulfanilic acid in the presence of a detrjent to from a red complex which is measured colorimetrically at 540 nm.

Specimen: Freshly prepared and nonhemolyzed serum.

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### Bilirubin (Direct Automated)

Method: colorimetric , End Point , Automated and manual.

Principle: Congugated (Direct) bilirubin react with diazotized sulfanilic acid to form a red complex in alkaline solution. The intensity of the color measured at 540 nm, is proportional to the concentration of direct bilirubin.

Specimen: Freshly prepared and nonhemolyzed serum.

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### Bilirubin (Manual)

Method: Colorimetric, End-point, Manual

Principle:

Conjugated (direct) bilirubin reacts directly with diazotized sulfanilic acid to form a red complex in alkaline solution. Direct bilirubin is, therefore, measured alorimetrically as a red azo dye at 546nm. The intensity of the color is proportional to the concentration of direct bilirubin. Unconjugated (indirect) bilirubin in the blood is protein-bound. In order to measure unconjugated bilirubin, therefore, it must be separated from the protein and made water-soluble by addition of the so-called "accelerator". In the presence of accelerator, conjugated and unconjugated bilirubin (i.e. total bilirubin) react with diazotized sulfanilic acid to form a green-colored azo-bilirubin complex. This reaction is carried out in the alkaline medium of Fehling Solution. The intensity of the green color is measured at 578nm and is proportional to the concentration of total bilirubin. Unconjugated (indirect) bilirubin is then calculated by difference.

Specimens:

Both serum and Plasma can be used. Samples should be protected from light as exposure to sunlight reduces the concentration of bilirubin. Although grossly hemolyzed samples should not be used, slight hemolysis may be acceptable.

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### Calcium

Method: Colorimetric, End-point , Automated/Manual

Calcium in the sample is permitted to form a colored complex with Arsenazo III in neutral pH. The intensity of the color is directly proportional to the concentration of calcium and is measured colorimetrically.

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### Calcium

Method: Colorimetric, End-Point, Manual

Principle:

Calcium in the sample is permitted to form a magenta-colored complex with cresolphthalein complexone (CPC) in alkaline solution. The intensity of the color is directly proportional to the concentration of calcium and is measured colorimetrically at 575nm.

Specimens:

Serum is the sample of choice. Plasma can only be used if it is heparinized since many other anticoagulants interfere with calcium ions. Hemolyzed samples should not be used. Urine samples, if not used fresh, should be acidified by treatment with concentrated hydrochloric acid in order to prevent precipitation of calcium salts.

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## Cholesterol

Method: Colorimetric, End-point, Manual/Automated

Principle:

The method used is based on the hydrolysis of cholesterol ester using the enzyme cholesterol esterase, to form nonesterified (free) cholesterol and free fatty acids. The non-esterified cholesterol is then oxidized, using cholesterol oxidase, to form hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) and cholest-4-en-3-one. In the presence of peroxidase the H<sub>2</sub>O<sub>2</sub> reacts with 4-chlorophenol and 4-aminophenazone to form a colored complex. The intensity of the color is directly proportional to the total cholesterol concentration in the sample and is measured colorimetrically at 520nm.

Specimens:

Either serum or plasma can be used. Presence of lipemia and gross hemolysis may interfere with results.

Contents:

Reagent 1: Enzyme Reagents Solution	3 x 100ml
Reagent 2: Cholesterol Standard (200mg/dl)	1 x 5ml
Number of Assays: Manual	300
Automated Dependent on Instrument Type	
Catalogue Number:	DFK-122

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Method: Colorimetric , Automated & Manual

Principle:

Copper in acidic solution is released from ceruloplasmin and reacts with 3,5-DiBr-PAESA to form a color complex. The intensity of the complex is proportional to copper concentration and is measured colorimetrically.

Specimen: Freshly prepared non-hemolyzed serum.

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## Creatinine

Method: Jaffe ,Two-point, Automated/Manual

Principle: Creatinine in the sample reacts with alkaline picrate to form a red colored complex. The rate of production of the complex is directly proportional to the concentration of creatinine. Specimens:

Serum or plasma can be used. Hemolyzed samples should be avoided. 24-urine creatinine can also be analyzed using this kit.

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## Creatinine

Method: Colorimetric, Two-point, Manual

Principle:

The method is based on the Jaffe reaction. Under alkaline conditions creatinine reacts directly with picrate ions to form a red-colored complex. The absorbance of the complex can be measured colorimetrically at 500nm and is directly proportional to the concentration of creatinine. Since the Jaffe reaction is not specific, the preliminary separation of creatinine from other serum components is necessary if the true creatinine level is to be determined. This is achieved by the use of an acid and a precipitating agent.

Specimens:

Serum or plasma can be used. Hemolyzed samples should be avoided. 24-hour urine creatinine can also be analyzed using this kit.

Contents:

Reagent 1: Precipitating Reagent	1 x 50 ml
Reagent 2: Acidic Reagent	1 x 100 ml
Reagent 3: Picric Acid	1 x 50 ml
Reagent 4: Alkaline Reagent	1 x 10 ml
Reagent 5: Creatinine Standard (10mg/dl)	1 x 10 ml
Number of Manual Assays:	100
Catalogue Number:	DFK-108

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### Creatinine Phosphokinase (CK)

Method: Enzymtic, UV Kinetic, Manual/Automated

Principle:

Measurement of CK is based on its ability to catalyze the reaction between creatine phosphate and adenine diphosphate (ADP) in acid solution. The adenosine triphosphate (ATP) thus formed then reacts with glucose to form glucose-6-phosphate (G-6-P). G-6-P then reacts with nicotinamide-adenine dinucleotide phosphate (NADP) and converts it to the reduced form (i.e. NADPH). The concentration of NADPH is directly proportional to the activity of CK and is monitored colorimetrically at 340nm. N-acetylcysteine is incorporated in to the reagents to activate CK. Diadenosine pentaphosphate and adenosine monophosphate (AMP) are incorporated to inhibit adenylate kinase (AK) which could interfere with the results.

Specimens:

Serum or heparinized plasma may be used. Grossly hemolyzed samples should avoided.

Contents:

Reagent 1: Buffer Solution	1 x 55ml
Reagent 2: Enzymatic Reagent	1 x 10 Vials
Number of Assays: Manual	50
Automated: Dependent on instrument type	
Catalogue Number:	DFK-145

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### Glucose

Method: Colorimetric, End-point, Enzymatic, Manual/Automated

Principle:

In the presence of oxygen, D-glucose is specifically converted to gluconic acid by the enzyme glucose oxidase (GOD). Other sugars do not take part in this reaction. Hydrogen peroxide, the other product of the reaction, is then mixed with a chromogen (i.e. 4-aminophenazone) in the presence of phenol and peroxidase to form a colored complex. The intensity of the color is directly proportional to the concentration of glucose and is determined colorimetrically at 500nm.

Specimens:

If immediate analysis is not possible, serum should be separated from cells immediately after drawing the blood since 7% of glucose content is lost per hour because of glycolysis caused by erythrocyte's enzymes. Hemolytic samples should be avoided. Plasma can also be used provided it is separated from the cells. Alternatively, sodium fluoride which acts as a preservative as well as an inhibitor of glycolytic enzymes can be added to the sample.

Contents:

Reagent 1: Enzymatic Reagent	3 x 100 ml
Reagent 2: Glucose Standard	1 x 7 ml
Number of Assays: Manual	300
Automated Dependent on Instrument type	
Catalogue Number:	DFK-120

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### Glucose 6-Phosphate Dehydrogenase, G6PD

Method: Qualitative, Manual

Principle:

G6PD catalyzes oxidation of gulcose 6-phosphate by NADP + produces NADPH and 6-phosphogluconolactone. NADPH reduces a redox dye, 2,6-dichlorophenolindophenol in the presence of phenazine methosulfate as an intermediate carrier. The rate of decolorization is proportional to the enzyme activity.

Specimens:

Whole blood collected on EDTA or heparin.

Contents:

Reagent 1: Buffer Solution	10ml
Reagent 2: Substrates	20 Vials
Reagent 3: Oxidizing Reagent	5ml
Reagent 4: Mineral Oil	15ml
Number of Assays: Manual	20
Catalogue Number:	DFK-117

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+Phosphogluconolactone + NADPH+H-۶ +Phosphate +NADP-۶ Glucose

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### Hemoglobin A1c Direct

HbA1c DIRECT (Immunoturbidimetry )

Method: Turbidimetry , Authomated and Manual

Principle:

This method is based on the interaction of antigen and antibody for direct determination of HbA1c in whole blood. Total hemoglobin and HbA1c have the same unspecific absorption rate to latex particles. Addition of mouse antihuman HbA1c monoclonal antibody forms latex-HbA1c-mouse antihuman HbA1c complex. When goat anti-mouse IgG polyclonal antibody is added to this complex agglutination is formed. The amount of the agglutination formed is proportional to the concentration of HbA1c absorbed on the surface of latex particles. The amount of the agglutination is measured colorimetrically.

Specimen: Whole blood samples collected with EDTA anticoagulant can be used. Prior to testing the sample should be mixed by gentle inversion to re-suspend settled erythrocytes.

Comtents:

Reagent R1:	Latex particles, Buffer, stabilizers.
Reagent R2:	Mouse anti-human HbA1c monoclonal antibody, gout antimouse IgG polyclonal antibody
Reagent R3 :	Lyse reagent, buffer, stabilizers.
Reagents volume: DFK-158:	28 ml; DFK- 159: 40 ml; DFK-160: 60 ml

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### High Density Lipoprotein (HDL) -Cholesterol

Method: Enzymatic, Direct Method, Manual/Automated

Principle: This method is a homogenous assay for direct measuring serum HDL-C level without the need for any pretreatment steps. First, substances with high affinity to LDL, VLDL and chylomicrons block them preventing not to participate in the enzymatic reaction. Then special surfactant that accelerates reagents of enzyme reaction with HDL-cholesterol is added.

Specimen: Freshly prepared serum or plasma (on EDTA) free of hemolysis can be used.

Reagents:

Reagent1(R1): Buffer pH 6.5 containing, blocking agents, surfactant, detergent, stabilizer.	1x45 ml
Reagent2(R2): Enzymatic reagents containing cholesterol esterase, cholesterol oxidase, peroxidase, detergents and stabilizers.	1x15 ml
Catalogue Number:	DFK-146

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### High Density Lipoprotein (HDL) Cholesterol

Method: Enzymatic, Colorimetric, End-point, Manual

Principle:

The measurement of HDL-cholesterol is based on precipitation. Phosphotungstic acid and magnesium ions, when added to the sample, precipitate the very-low-density-lipoprotein (VLDL) and the low-density-lipoprotein (LDL) leaving only the HDL-cholesterol in the supernatant. The supernatant is then analyzed colorimetrically (520nm) using total cholesterol reagents. The measured cholesterol concentration refers to HDL-cholesterol.

Specimens:

Either serum or plasma can be used. Presence of lipemia and gross hemolysis interferes with results.

Contents:

Reagent 1: Precipitating Solution	3 x 50ml
Reagent 2: Cholesterol Standard (20 mg/dl)	1 x 5ml
Number of Assays: Manual	300
Catalogue Number:	DFK-130

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### LDL-Cholesterol , Direct

Method: Enzymatic, Automated and manual

Principle:

In this method non-LDL lipoprotein particles (CM,HDL,VLDL) are released by a specific detergent and their cholesterol will be used up by enzymatic reagents and be in a non-color forming reaction without chromogen coupler. In the next step addition of another detergent solubilized LDL particle and the released cholesterol and cholesterol esters oxidized by cholesterol oxidase and react with chromogenic coupler for the color formation . The intensity of the color is proportional to LDL concentration.

Specimens: Freshly prepared serum or plasma (on EDTA or heparin) may be used. Serum or plasma should be separated from coagulum as soon as possible (within 2 hour) and not to be stored for more than 12 hours at room temperature or 5 days at refrigerator.

Contents:

Reagent R1: Cholesterol esterase, cholesterol oxidase, per oxidase, detergent, chromogen, stabilizers.	45 ml
Reagent R2: Detergent, chromogen, stabilizers.	15 ml
Number of assays:	Dependent on the instrument parameters.
Catalogue Number:	DFK-157

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### Serum Iron Ferrozine Method

Method: Colorimetric, End-point, Manual/Automated

Principle:

Transferrin – bound iron is released into an acidic buffer solution (pH 4.5) without precipitation of proteins. The released iron is reduced and subsequently reacts with Ferrozine to form a water-soluble magenta complex. The intensity of the color is directly proportional to the concentration of iron and is measured colorimetrically at 562nm.

Specimens:

Serum and sample of choice. Plasma should not be used because of the precipitation reaction of iron with most anticoagulants. Serum should be free of hemolysis.

Contents:

Reagent 1: Buffer Solution	2 x 50ml
Reagent 2: Chromogenic Reagent	2 x 10ml
Reagent 3: Iron Standard (200 µg/dl)	1 x 10ml
Number of Assays: Manual	100
Automated Dependent on Instrument type	
Catalogue Number:	DFK-121

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### TIBC Direct Colorimetric

Method : Automated

Principle : Acidic solution of iron is first added to serum then iron plus the iron released from transferrin and that of serum will all bind to the color Reagent. Then the pH at the solution will be brought to neutral where the iron will be bind to transferrin and the absorbance at 660nm will be decreased which is proportional to TIBC Concentration .

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Method: Manual

Principle:

Iron in serum is bound to transferrin. Normally, however, only about one-third of the total saturation capacity of transferrin is bound to iron. In order to determine the total iron-binding capacity of transferrin (TIBC), ferric iron is added in excess to the sample, thus saturating all the transferrin with iron. Any unbound ferric iron is then removed by the addition of magnesium carbonate. TIBC is then determined calorimetrically by using a serum iron kit. If a serum iron determination is performed at the same time, the difference between the TIBC value and serum iron represents the transferrin that was not bound to iron. This quantity is known as unsaturated iron-binding capacity (UIBC).

Specimens:

Serum is the sample of choice. Plasma should not be used because of the precipitation reaction of iron with most anticoagulants. Serum should be free of hemolysis.

Contents:

Reagent 1: Iron Solution	1 x 55ml
Reagent 2: Magnesium Carbonate(Powder)	1 X 3g
Number of Manual Assays:	50
Catalogue Number:	DFK-102

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### Lactate Dehydrogenase (LDH)

Method:  
Colorimetric, UV Kinetic, Manual/Automated  
Principle:

All five LDH isoenzymes catalyze the following reversible reaction:

Pyruvate + NADH → Lactate + NAD<sup>+</sup>

–The reaction equilibrium at pH 7.5 favors the reduction of pyruvate to lactate at a rate that is directly proportional to the activity of LDH. NADH has a concentration – dependent absorbance at 340nm whereas the other reagents do not. The reductions in absorption of NADH is proportional to the activity of LDH and is measured colorimetrically. Total LDH activity is determined in this method.

Specimens:

Serum and sample of choice. Plasma can also be used provided that heparin, citrate and oxalate are not used as anticoagulants. Hemolyzed samples must be avoided.

Contents:

Reagent 1: Pyruvate- Buffer Reagent	2 x 40 ml
Reagent 2: Co-enzyme (NADH)	1 x 20 ml
Number of Assays: Manual	100
Automated Dependent on Instrument type	
Catalogue Number:	DFK-131

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### Magnesium

Method: Colorimetric, Automated & Manual  
Principle:

Magnesium ions in alkaline solution react with Xylidyle Blue to form a color complex. The intensity of the color is proportional to magnesium concentration and is assayed by colorimetric measurement.

Specimens:

- 1- Serum or plasma (on heparin) . Hemolysis, EDTA, oxalate and Sodium fluoride will affect the test results.
- 2- 24 h Urine . Urine pH should be adjusted to 3-4 with HCl and diluted 1:5 with distilled water prior to use.

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Method: UV method, Automated/Manual

Principle:

Inorganic phosphate can be quantified by mixing the sample with molybdate ions (in acid solution) to form a phosphomolybdate complex. This complex can be analyzed at 340 nm. The absorbance of the complex is proportional to phosphorus concentration.

Specimens:

Either serum or plasma prepared on heparin can be used. The concentration of phosphate, however, is slightly higher in plasma as compared with serum. Even slightly hemolyzed samples should not be used. 24-hour urine samples may also be analyzed by this method.

Contents:

Reagent 1: Ammonium molybdate reagent	2 x 50 ml
Reagent 2: Standard phosphorus (5 mg/dl)	1 x 5 ml
Number of tests: Manual : 100 ; Automated : Varies with instrument type.	
Catalogue Number:	162

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## Phosphorus

Method: Colorimetric, End-point, Manual

Principle:

Phosphorus in the blood occurs as the phosphate, either free (inorganic phosphate) or in combination with an organic group (organic phosphate). The latter is mainly confined to erythrocytes. Therefore, phosphate determinations in serum refer mainly to inorganic phosphates. Inorganic phosphates can be quantified by mixing the sample with molybdate ions (in acid solution) to form a phosphomolybdate. Addition of a mild reducing agent will then reduce the phosphomolybdate to molybdenum blue, which can be analyzed colorimetrically at 640nm. The intensity of the color is directly proportional to the concentration of inorganic phosphates.

Specimens:

Either serum or plasma can be used; however, the concentration of phosphate is slightly higher in plasma as compared with serum. Serum/plasma should be separated from cells as soon as possible. Even slightly hemolyzed samples should not be used.

Contents:

Reagent 1: Reducing Agent	1 x 110 ml
Reagent 2: Molybdate ion Reagent	1 x 25 ml
Reagent 3: Protein Solvent	2 x 80 ml
Reagent 4: Phosphorus Standard (5mg/dl)	1 x 10 ml
Number of Assays:	100
Catalogue Number:	DFK-110

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### Total Protein

Method: Colorimetric, End-point, Manual/Automated

Principle:

The biuret method is used to determine total protein levels. The reaction between proteins and cupric ions forms a violet-colored complex. The intensity of the color is measured colorimetrically (546nm) and is directly proportional to the concentration of proteins.

Specimens:

Serum in the sample of choice since the fibrinogen present in plasma is a protein and interferes with the analysis. Hemolyzed serum should not be used.

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### Triglycerides GPO-PAP Method

Method: Enzymatic, Colorimetric, End-point, Manual/Automated

Principle:

Triglycerides are hydrolysed to glycerol and free fatty acids by the action of lipase. The liberated glycerol is then phosphorylated by adenosine triphosphate (ATP) in the presence of glycerol kinase to form glycerol-3-phosphate (G-3-P) and adenine diphosphate. G-3-P is then oxidized by glycerol-3-phosphate oxidase to form dihydroxy acetone phosphate (DHAP) and hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>). H<sub>2</sub>O<sub>2</sub> is then mixed with a chromogen (4-aminophenazone) in the presence of p-chlorophenol to form a colored complex. The intensity of the color is measured colorimetrically at 505nm and is directly proportional to the concentration of triglycerides.

Specimens:

Serum or plasma may be used. Hemolyzed samples should not be used.

Contents:

Reagent 1: Enzyme solution	3 x 100 ml
Reagent 2: Standard (200 mg/dl)	1 x 10 ml
Number of Assays: Manual	300
Automated Dependent on Instrument Type	
Catalogue Number:	DFK-124

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### UREA Berthelot Method , Enzymatic (Automated&Manual)

Method: Berthelot Method, Manual

Principle:

Urea is converted to ammonium ion and carbondioxide by urease. The ammonium ion produces a complex with salicylate and hypochlorite ions which absorbs the light at 578 nm. The absorbance is proportional to the urea concentration.

Specimens:

Either serum or plasma can be used for analysis. If plasma is to be used anticoagulants containing ammonium or fluoride ions must be avoided. This kit can also be used to determine the concentration of urea in urine.

Contents:

Reagent 1: Buffer Solution (pH 7.2)	2 x 50 ml
Reagent 2: Enzyme Solution	1 x 2 ml
Reagent 3: Sodium-hypochlorite	2 x 50 ml
Reagent 4: Standard Solution (50 mg/dl)	5 ml
Number of Assays:	200
Catalogue Number:	DFK-116

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### UREA Urease-UV Method (Automated&Manual)

Method: Enzymatic, UV, Two-point, Manual/ Automated

Principle:

Urea is hydrolyzed to ammonium ions by urease. These ions then react with α-ketoglutarate in the presence of glutamate dehydrogenase and nicotinamide adenine diphosphate (NADH). This reaction results in the oxidation of NADH to NAD<sup>+</sup>, which in turn decreases the spectrophotometric absorption (340nm) of the solution. The decrease in absorption is proportional to the concentration of urea.

Specimens:

Either serum or plasma can be used for analysis. However, if plasma is to be used anticoagulants containing ammonium or fluoride ions must be avoided. This kit can also be used to determine the concentration of urea in urine.

Contents:

Reagent 1: Tris Buffer (pH 7.8)	3 x 80 ml
Reagent 2: Substrate Reagents	1 x 60 ml
Reagent 3: Urea Standard (50 mg/dl)	1 x 10 ml
Number of Assays: Manual	300
Automated Dependent on Instrument Type	
Catalogue Number:	DFK-127

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Method: Colorimetric, End-Point, Manual

Principle:

Diacetyl monoxime reacts with water in an acidic solution to form Diacetyl. Urea reacts with diacetyl in the presence of thiosemicarbazide to form a pink-colored diazine derivative. The intensity of the color is directly proportional to the concentration of BUN and is monitored calorimetrically at 520nm.

Specimens:

It is advisable to use serum, however plasma containing certain anticoagulants (e.g. EDTA or heparin) may also be used.

Contents:

Reagent 1: Chromogenic Reagent	2 x 100ml
Reagent 2: Acidic Reagent	2 x 100ml
Reagent 3: BUN Standard (16 mg/dl)	1 x 5ml
Number of Manual Assays:	400
Catalogue Number:	DFK-104

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### Uric Acid Uricase / PAP

Method: Enzymatic Colorimetric, Uricase/PAP, Manual/ Automated

Principle:

The method is based the oxidation of uric acid to allantoin and hydrogen peroxide by uricase. Hydrogen peroxide in the presence of peroxidase reacts with 4-aminophenazone and 3,5, dichloro benzene sulfonate (DHBS) to form a quinoneimine dye. The concentration of the dye is directly proportional to the uric acid concentration which measured colorimetrically.

Specimens:

Either serum, plasma prepared on heparine or a 1:10 dilution urine samples may be used.

Contents:

Reagent 1: Phosphate Buffer pH 7.8	2 x 100 ml
Reagent 2: Enzyme Reagents	2 x 5 ml
Reagent 3: Standard Solution (6 mg/dl)	1 x 7 ml
Number of Assays: Manual	200
Automated: Varies with the type of instrument	
Catalogue Number:	DFK-133

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### Uric Acid (Manual)

Method: Colorimetric, End-point, Manual

Principle:

Uric acid is oxidized to allantoin by a phosphotungstic acid reagent in alkaline solution. Phosphotungstic acid is reduced in this reaction to tungsten blue, which is monitored calorimetrically at 640nm. The intensity of the color is directly proportional to the concentration of uric acid.

Specimens:

Either serum or plasma can be used as long as they are free from hemolysis. This kit can also be used to determine the concentration of uric acid in urine.

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### Urine and CSF Proteins

Methods: Micro - Turbidimetric / Automated

Principle:

Proteins present in the urine or CSF react with denaturing reagents in alkaline solution to produce a turbid solution. The extent of turbidity is proportional to the protein concentration and is measured colorimetrically at 546 nm.

Specimens:

1. Twenty four – hour urine sample collected in a bottle containine 10 ml concentrated HCL.
2. Freshly prepared CSF

Contents:

Reagent 1 : Alkaline solution	2x45 ml
Reagent 2 : Denaturing reagent	2x15 ml
Reagent 3 : Protein standatd (200mg/dl)	1x5 ml
Automated: Dependent on the instrument type.	
Catalogue Number:	DFK-155

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### Urine and CSF proteins (Pyrogallol Red)

Method: Colorimetric, End -point, Automated I manual

Principles:

Proteins present in the urine or CSF reacts with pyrogallol red and molybdate in acidic solution to form a colored - complex. The intensity of the color is directly proportional to the concentration of protein and is measured colorimetrically at 600 nm.

Specimens:

1. Twenty four- hour urine sample collected in a bottle containing 10 ml concentrated HCL.
2. Freshly prepared CSF

Contents:

Reagent 1: Pyrogallol red solution	3 x 50 ml
Reagent 2: Standard protein (30 mg/dl)	1 x 5 ml
Number of Assays:	
Manual:	150
Automated: Dependent on the instrument type.	
Catalogue Number:	DFK-106

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### Urinary Citrate

Method : Enzymatic , Automated & manual

Principle : Citrate is Converted to Oxaloacetate and acetate by citrate lyase . Oxaloacetate then reacts with NADH by Malate Dehydrogenase producing malate and NAD<sup>+</sup> .The Change in NADH concentration is proportional to the citrate which is measured at 340 nm.

Specimen : Twenty four-hour Urine sample collected in a bottle containing Concentrated HCL can be used.

Contents :

Reagent 1: Tris – HCL Buffer	9 ml
Reagent 2: Substrate	2.5 ml
Reagent 3: Lyophilized Enzyme	5 Vials
Reagent 4: Citrate Standard (200mg/L)	5 ml
Number of Assays: Dependent of the instrument type.	
Catalogue Number:	151

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Method: Enzymatic, Manual

Principles:

Citrate Lyase catalyzes conversion of citrate to oxaloacetate and acetate. Same oxaloacetate may be decarboxylase which possibly contaminates the citrate lyase preparation. Both  $\alpha$ -ketoacids react with phenylhydrazine at acidic pH and produce phenylhydrazone which absorbs the light at 340 nm. The absorbance is directly proportional to citrate concentration.

Specimens:

Twenty four- hour urine sample collected in a bottle containing 10 ml concentrated HCL can be used.

Contents:

Reagent 1: Tris- HCl Buffer	25 ml
Reagent 2: Lyophilized Enzyme	5 Vials
Reagent 3: Phenylhydrazine	1 ml
Reagent 4: Citrate Standard (200 mM)	5 ml
Number of Assays: Manual	20
Catalogue Number:	DFK-105

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### Urinary Oxalate Enzymatic Method (Automated & Manual)

Method: Manual/Automated, Enzymatic

Principles:

Urine oxalate is converted to carbon dioxide and hydrogen peroxide. Hydrogen peroxide in the presence of peroxidase oxidizes a chromogen to form a colored complex. The intensity of the color is directly proportional to the concentration of oxalate and is determined calorimetrically at 590 nm.

Specimens:

Twenty four- hour urine sample collected in a bottle containing 10 ml concentrated HCL can be used.

Contents:

Reagent 1: Buffer Solution	25 ml
Reagent 2: Lyophilized Enzyme	5 Vials
Reagent 3: Active Charcoal	12 g
Reagent 4: Oxalate Standard (0.25 mM)	5 ml
Reagent 5: Sample Diluent	60 ml
Reagent 6: Lyophilized Color Reagent	1 Vials
Number of Assays: Manual	20
Automated:	60
Catalogue Number:	DFK-103

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## Urinary Calculi (Stone analysis)

### Calculus Tests Contained In This Kit:

A complete ready-for-use reagent kit for semi-quantitative determination of the most important components of urinary calculi, i.e. uric acid, cystine, phosphate, agnesium, ammonium, calcium and oxalate.

Method: Chemical, Manual

### Principles:

Uric acid is detected by reduction of alkaline phosphotungstate to tungsten blue; cystine by the red color formed in the cyanide-nitroprusside test; phosphate by the blue color formed by reduction of molybdato-phosphoric acid to molybdenum blue; magnesium by the red color developed with xylidyl blue; ammonium ions are assayed using Nessler's reaction. Calcium is measured by titration against Titriplex III; oxalate by the color formed in reacting with iron chloride and sulphosalicylic acid. The various salts of particular ions can then be determined using the specially-calibrated ruler provided.

### Specimens:

The excreted stone must be crushed to a homogeneous powder, which is then made in to a solution and used in all subsequent tests.

### Contents:

Reagent 1: Sulphuric Acid Solution	2 x 11 ml
Reagent 2: Sodium Hydroxide Solution	1 x 17 ml
Reagent 3: Calcon Carboxylic Acid Solution	1 X 3.5 g
Reagent 4: Titriplex III Solution	1 x 18 ml
Reagent 5: Borate Buffer Solution	1 x 18 ml
Reagent 6: Iron Chloride Solution	1 x 12 ml
Reagent 7: Sulphosalicylic Acid Solution	1 x 12 ml
Reagent 8: Potassium Mercury Iodide Solution	1 x 12 ml
Reagent 9: Ammonium Molybdate	1 x 20 ml
Reagent 10: Reduction Solution	1 x 20 ml
Reagent 11: Borate Buffer Solution	1 x 20 ml
Reagent 12: Color Reagent	1 x 20 ml
Reagent 13: Molydatophosphoric Acid	1 x 12 ml
Reagent 14: Ammonia Solution	1 x 20 ml
Reagent 15: Sodium Sulphite	1 x40 g
Reagent 16: Sodium Nitroprusside	1 x40 g
A control sample is also enclosed	
Number of Manual Assays:	80
Catalogue Number:	DFK-107

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## Urinary Cystine

Method : Colorimetric , End Point , Manual

Principle : Cystine in urine is reduced to cysteine and cysteine sulfonate. The reduced substances then react with phosphotungstic acid in acidic solution to form a blue color complex. The intensity of the color is directly proportional to the Cystine concentration and is measured at 600 nm.

Specimens : Twenty four-hour urin sample may be used.

### Contents :

Reagent 1 : Acetate buffer	1 × 50 ml
Reagent 2: Reducing reagent	1 × 50 ml
Reagent 3 : Mercury chloride solution	1 × 10 ml
Reagent 4: Colorimetric reagent	1 × 50 ml
Reagent 5: Cystine standard (200 mg/l)	1 × 5 ml
Number of Assays :	
Manual	60 tests
Automated:	Dependent on instrument type
Catalogue number:	DFK-152

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## ZINC

Method: Colorimetric, End Point, Automated and Manual

Principle: Ionic form of zinc react with 5-Br-PAPS to form a red color complex. The intensity of the color is directly proportional to zinc concentration and is measured at 560 nm.

Specimen: Non-hemolyzed serum or plasma can be used. Twenty-four urine sample may be also used.

### Reagents:

Reagents 1: Chromogenic reagent	1 × 50 ml
Reagents 2: Zinc standard (200 µg/dl)	1 × 5 ml
Number of Assays:	
Manual:	50 tests
Automated:	Dependent on instrument type.
Catalogue number:	DFK-154

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Hematology Solution

General Diluent	MS9, ERMA,HYCEL.....
Sysmex Diluent	Kx21, K-1000 , K-800
Lysis	Kx-21
Lysis of WBC	K-1000
Lysis of Hb	K-1000
Enzymatic Cleanser	ERMA , HYCEL.....

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